

Certificate of Conformity

Konformitätsbescheinigung
Certificat de Conformité
Certificato di conformità
Certificado de conformidad
Certificado de Conformidade



Page 1 of 1

Part Number	8080A-60	8080L-60
	8080A-150	8080L-150
	8080A-450	8080L-450
	8080A-900	8080L-900
	8080A-1200	8080L-1200
	8080A-5500	
	8080A-10500	
Description	Sterile Bags	
Material	LDPE (Low Density Polyethylene)	
Temperature Rating	80°C to – 100°C	
Free from:	Pyrogens, R-Nase & D-Nase	
Shelf Life:	5 years from date of manufacture	
Date	21 January 2025	

We confirm that these bags are manufactured from LDPE which conforms to FDA 21 CFR 177.1520.

The manufacturer confirms these bags are manufactured from fully synthetic polymers. No animal sourced materials are used in the production of PureBags.

The bags are made from extruded virgin polyethylene at temperatures up to 200°C. The polyethylene tube manufacturing process is done by extrusion and the internal tubing is never exposed to the external environment. This ensures the bags sterility up until the point at which it is opened by the user.

The sterility of the manufactured bags is controlled by periodically sampling and analysing various production lines. In addition, sterilization by the extrusion process has been scientifically validated. Independent testing has confirmed that the sterile sampling bags we supply remain sterile for at least five years if they are not opened.

The information detailed above has been given to us by our suppliers, it is passed on in good faith and is to the best of our knowledge correct.

For and on behalf of Sampling Systems Ltd.

A blue ink signature of Kuldeep Flora, written in a cursive style.

Kuldeep Flora
Quality Manager

This certificate is valid for 3 years from the date given above. Download the most recent certificate at: www.sampling.com/certificates



Sampling Systems Ltd.

2 Trillennium Point, Gorsey Lane, Coleshill. B46 1JU. UK.
Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com