

1. Part Number: 1035A-1000

2. Photo:



3. Product Information:

Description	Sample Sleeve Sampler Handle 1000mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Overall Length:	1070mm
Weight:	850g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Nov-24