

**1. Part Number:** 1035A-1500

**2. Photo:**



**3. Product Information:**

Description	Sample Sleeve Sampler Handle 1500mm
Materials:	316L stainless steel
External Polish:	Better than 0.8 microns Ra
Overall Length:	1580mm
Weight:	1130g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Nov-24