

1. Part Number: 1070A-350

2. Photograph:



3. Product Information:

Description	MicroMate 350mm
Materials of Construction:	316L stainless steel Nylon Clutch (non product contact part)
Nominal Sample Volume:	Variable up to 0.1ml
Overall Length:	500mm
Length of Sampling Tube:	350mm
Diameter of Sampling Tube:	5mm
Weight:	200g
Finish:	<0.8 microns Ra
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 12-Nov-24