

1. Part Number: 1070A-500

2. Photograph:



3. Product Information:

Description	MicroMate 500mm
Materials of Construction:	316L stainless steel Nylon Clutch (non product contact part)
Nominal Sample Volume:	Variable up to 0.1ml
Overall Length:	650mm
Length of Sampling Tube:	500mm
Diameter of Sampling Tube:	5mm
Weight:	220g
Finish:	<0.8 microns Ra
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are free of animal products
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 12-Nov-24