

1. Part Number: 1290K-1000


2. Photograph:



Note- 3 sample inserts are shown for information only – they must be ordered separately

3. Product Information:

Description:	Multi Thief 1000mm
Method of Construction:	Fully welded handles – screw in tip
Finish	External Bright Polished Internal Machine finish
Nominal Length:	1060mm
Diameter of Outer Tube:	40mm
No. of Sampling Slots	3 off
Position of Sampling Slots:	



No. of inserts required	3 off
Nominal Weight of Sampler:	3.4 Kg (sampler only – no inserts)
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Nov-24