

1. Part Number: 1320A-01

2. Photo:



3. Product Information:

Description	Sack Master – with Plain End
Method of Construction:	Single piece construction
Materials:	316L stainless steel
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	400mm
Outer Tube Diameter:	25mm
Inner Diameter:	22mm
Nominal Weight:	335g
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 14-Nov-24