

- 1. Part Number:** 2000A-1000
2000T-025
2000T-050
2000T-075
2000T-100

2. Photo:



3. Product Information:

Description: U-D Cohesive Sampler & Tips
Materials: 316L stainless steel
(Note: the stainless steel spring is not grade 316L)
Method of Construction: Crevice free, all welds ground and polished

Sampler Body (Part No: 2000A-1000)

Max. Diameter: 26mm
Overall Length (sampler & tip): 1080mm
External Finish: <1 Micron Ra
Nominal Weight of Sampler (2000A-1000): 660g

Tips

Part No.	Length	Overall Diameter	Depth of Sample Chamber	Weight
2000T-025	30mm	26mm	3.2mm	60g
2000T-050	30mm	26mm	5.5mm	57g
2000T-075	30mm	26mm	7.6mm	56g
2000T-100	35mm	26mm	9.8mm	74g

BSE/TSE statement: All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 19-Nov-24