


**1. Part Number:** 2015A-25

**2. Photo:**



### 3. Product Information:

Description	Cohesive Sampler 25ml
Materials:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.8 microns Ra
Nominal Sample Volume:	25ml
Overall Length:	1030mm
Diameter of Outer Tube:	25mm
Scraper Length	110mm
Nominal Weight of Sampler:	1.9Kg
BSE/TSE statement:	All polishing compounds used in the manufacture of this sampler are of vegetable origin and so are BSE/TSE free.
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 19-Nov-24