

1. Part Number: 4000A-01

2. Photograph:



3. Product Information:

Description:	Multilevel Tablet and Capsule Sampler (removeable cups)
Material of Construction:	Frame: 316L stainless Steel Cups: 304 Stainless Steel
Method of Construction:	Cups lift out of frame
Finish	Ra <1.0µm
Number of Sampling Cups:	3
Volume of Sampling Cups:	20ml
Overall Length:	800mm
Nominal Weight:	800g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 29-Nov-24