

**1. Part Number:** 5020T-01

**2. Photo:**



**3. Product Information:**

Description Pump Sampler Set – PTFE 180ml

Set Contents:

Part No.	Description	Quantity	Material
5020A-09	Hand Pump	1	Plastic with polymer O-ring seals
5020A-11	Bottle Adapter*	1	PTFE
5020T-05	Sampling Tube (OD 8mm) x 5m length)	1	PTFE
5020A-07	Tube Cutter	1	Plastic with metal blade
5020A-04	Sample Bottle 180ml	1	PFA
5020A-98	Tube Weight	1	Stainless Steel

\* The Hand Pump is suitable for use with tubes OD from 5mm to 8mm

\*Sampling Systems offer a PTFE bottle adapter that can be custom made to suit any bottle with a screw thread

Recommended Storage Conditions: Dry and ambient temperature

Recommended Spares:

5020T-05	PTFE Tubing – 5m length
5020T-100	PTFE Tubing – 100m length
5020T-100D	PTFE Tubing – 100m length & documentation**
5020T-03	Sample Bottle 90ml

\*\* Documentation Package includes a certification of FDA compliance, date of manufacture.

Notes: 1. Once the sample is full the operator should stop operating the pump. If the user continues to operate the pump then the liquid will be drawn into the body of the pump and will leak out.  
2. This pump must not be used with flammable solvents

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Nov-24