

1. Part Number: 5021P-250

2. Photo:



3. Product Information:

Description Pump Sampler Set – PE 250ml

Set Contents:

Part No.	Description	Quantity	Material
5020A-09	Hand Pump	1	Plastic with polymer O-ring seals
5020P-10	Sampling Tube (OD 6mm, ID 4.5mm, Length 13m)	1	Polyethylene
5020A-07	Tube Cutter	1	Plastic with metal blade
5021H-250	Sample Bottle 250ml	5	HDPE

* The Hand Pump is suitable for use with tubes OD from 5mm to 8mm

Recommended Storage Conditions: Dry and ambient temperature

Dimensions of Case: 390 x 300 x 100mm

Weight of Set: 1.1 kg

Notes: 1. Once the sample is full the operator should stop operating the pump. If the user continues to operate the pump then the liquid will be drawn into the body of the pump and will leak out.
2. This pump must not be used with flammable solvents

4. Approval:

Signed:



Name: Kuldeep Flora
 Position: Quality Manager
 Date: 29-Nov-24