


1. Part Numbers: 5045A-1000

2. Photograph:



3. Product Information:

Description Multi Dipper 2 – 1000mm

Material of Construction:


Frame:	316L stainless steel
Chain:	316 stainless steel
Shackles:	316 stainless steel
Bottle Lids:	316L stainless steel
Bottles (8165A-250):	HDPE

Max. Width: 115mm

Overall Length (sampler & bottles): 1020mm

No. of Sample Bottles: 3

Bottle Positions:



External Finish: <1 Micron Ra

Nominal Weight: 2.1 Kg (sampler & empty bottles)

Recommended Storage Conditions: Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 21-Jan-25