


1. Part Numbers: 5072P-500

2. Photograph:



3. Product Information:

Description	Sub Sampler Pro PTFE 500ml
Material of Construction:	PTFE
Method of Construction:	All parts screw together
Nominal Sample Volume:	500ml
Overall Length:	371mm
Maximum Width:	70mm
Nominal Weight (without product):	1070g
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Dec-24