


1. Part Number: 5100A-100

2. Photo:



3. Product Information:

Description	PTFE Dipper 100ml
Product Contact Parts:	PTFE (Polytetrafluoroethylene) <ul style="list-style-type: none"> • Pure, virgin PTFE • USP Class VI compliant • FDA conforming
Nominal Capacity:	100ml
Diameter of Cup:	55mm
Handle Length:	600mm
Overall Length:	690mm
Working Temperature Range:	-200 to +280°C
Nominal Weight:	435g
Recommended Storage Conditions:	Dry and ambient temperature

Notes: Increase the effective reach of this sampler by using Extension Rods

- PTFE Extension Rod 600mm Part No. 5100E-600
- PTFE Extension Rod 1000mm Part No. 5100E-1000

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 12-Dec-24