


**1. Part Number:** 5100A-500

**2. Photo:**



**3. Product Information:**

|                                 |  |
|---------------------------------|--|
| Description                     | PTFE Dipper 500ml  |
| Product Contact Parts:          | PTFE (Polytetrafluoroethylene) <ul style="list-style-type: none"><li>• Pure, virgin PTFE</li><li>• USP Class VI compliant</li><li>• FDA conforming</li></ul> |
| Nominal Capacity:               | 500ml  |
| Diameter of Cup:                | 80mm   |
| Handle Length:                  | 600mm  |
| Overall Length:                 | 740mm  |
| Working Temperature Range:      | -200 to +280°C   |
| Nominal Weight:                 | 655g   |
| Recommended Storage Conditions: | Dry and ambient temperature  |

Notes: Increase the effective reach of this sampler by using Extension Rods

- PTFE Extension Rod 600mm Part No. 5100E-600
- PTFE Extension Rod 1000mm Part No. 5100E-1000

**4. Approval:**

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 12-Dec-24