

1. Part Numbers: 5132B-300

2. Photograph:



3. Product Information:

| | |
|---------------------------------|----------------------------------|
| Description | Immersion Cylinder BRASS – 300ml |
| Material of Construction: | Brass |
| Nominal Volume: | 300ml |
| Body Length | 125mm |
| Body Diameter: | 70mm |
| Nominal Weight: | 1.2kg |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 12-Dec-24