


**1. Part Number:** 7000T-2000

**2. Photo:**



**3. Product Information:**

Description	Viscous Sampler PTFE/FEP 2000mm
Materials:	PTFE & FEP (for product contact parts)
Max Sample Volume:	650ml (approx.)
Diameter of Outer Tube:	25mm
Max. Diameter of Tip:	33mm
Overall Length:	2150mm
Nominal Weight:	1415g
Recommended Storage Conditions:	Dry and ambient temperature

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24