

1. Part Number: 8005A-01

2. Photograph:



3. Product Information:

Description:	Sealing Sleeve 25mm
Use:	For use with 25mm diameter Pocket Samplers
Material of Construction:	HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free
Method of Construction:	Single piece
Max. Sampler Diameter:	25mm
Sleeve Length:	150mm
Nominal Weight:	11g
Number of Sleeves per Pack:	3
Recommended Storage Conditions:	Dry and ambient temperature

Note: The Pocket Sampler shown in the photograph is for illustrative purposes only. The sampler must be purchased separately to the Sealing Sleeves

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-24