


1. Part Number: 8025H-115S

2. Photograph:



3. Product Information:

Description:	SteriWare® Tweezers - Sterile
Use:	Single use, disposable tweezers (forceps)
Material of Construction:	High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> • Conforms to FDA CFR 177.1520 • EU Regulations 10/2011 compliant • EC regulation 1935/2004 compliant • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Overall Length:	115mm
Nominal Weight:	3g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 14-Jan-25

IMPORTANT NOTE

This product is going to be replaced by [SteriWare Forceps](#)



Click on the links below to find out about the replacement products

Old Part Number	New Part Number
8025H-115	8630A-115
8025H-115S	8630A-115S