


**1. Part Number:** 8050H-500S

**2. Photograph:**



### 3. Product Information:

|                                         |                                                                                                                                                                                                                                                                                                       |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description:                            | SteriWare® LiquiThief - 500mm - Sterile                                                                                                                                                                                                                                                               |
| Use:                                    | Single use, disposable liquid sampler                                                                                                                                                                                                                                                                 |
| Material of Construction:               | High Density Polyethylene (HDPE) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to REACH regulation 1907/2006</li> <li>• BSE/TSE free</li> </ul> |
|                                         | O-ring Seal – virgin Silicone Rubber <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.2600</li> </ul>                                                                                                                                                                                 |
| Method of Construction:                 | Ultrasonic welding                                                                                                                                                                                                                                                                                    |
| Overall Length                          |                                                                                                                                                                                                                                                                                                       |
| Fully closed:                           | 524mm                                                                                                                                                                                                                                                                                                 |
| Fully extended:                         | 885mm                                                                                                                                                                                                                                                                                                 |
| Diameter of Outer Tube:                 | 21mm                                                                                                                                                                                                                                                                                                  |
| Diameter of hole in tip:                | 1.5mm                                                                                                                                                                                                                                                                                                 |
| Width of Handle:                        | 60mm                                                                                                                                                                                                                                                                                                  |
| Nominal Weight of Sampler:              | 82g                                                                                                                                                                                                                                                                                                   |
| Moulding Environment:                   | Class 100,000 Cleanroom                                                                                                                                                                                                                                                                               |
| Assembly & Primary Packing Environment: | Class 10,000 Cleanroom                                                                                                                                                                                                                                                                                |
| Individually Bagged?                    | Yes (heat sealed PE/PET bag)                                                                                                                                                                                                                                                                          |
| Method of Sterilisation:                | Gamma Irradiation (25 to 45 kGy)                                                                                                                                                                                                                                                                      |

Quantity per Box: 20

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 2-Apr-25