

**1. Part Number:** 8063H-500S

**2. Photograph:**



### 3. Product Information:

|                                        |                                                                                                                                                                                                                                                  |
|----------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description:                           | SteriWare® Powder Lance - 500mm - Sterile                                                                                                                                                                                                        |
| Use:                                   | Single use, disposable powder sampler                                                                                                                                                                                                            |
| Material of Construction:              | HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:                | Ultrasonic welding                                                                                                                                                                                                                               |
| Nominal Length:                        | 615mm                                                                                                                                                                                                                                            |
| Blade Length:                          | 500mm                                                                                                                                                                                                                                            |
| Diameter of Blade:                     | 22mm                                                                                                                                                                                                                                             |
| Nominal Weight of Sampler:             | 85g                                                                                                                                                                                                                                              |
| Moulding Environment:                  | Class 100,000 Cleanroom                                                                                                                                                                                                                          |
| Welding & Primary Packing Environment: | Class 10,000 Cleanroom                                                                                                                                                                                                                           |
| Individually Bagged?                   | Yes (heat sealed PE/PET bag)                                                                                                                                                                                                                     |
| Method of Sterilisation:               | Gamma Irradiation (25-45kGy)                                                                                                                                                                                                                     |
| Quantity per Box:                      | 20                                                                                                                                                                                                                                               |
| Recommended Storage Conditions:        | Dry and ambient temperature                                                                                                                                                                                                                      |
| Shelf Life:                            | 5 years from the date of manufacture                                                                                                                                                                                                             |
| ATEX:                                  | This device has not been ATEX tested                                                                                                                                                                                                             |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Jun-25