

1. Part Number: 8080L-1200**2. Product Information:**

Description: PureBag with Labelling Area 1200ml - Sterile

Use: Single use, disposable sampling bags

Material of Construction: Light Density Polyethylene (LDPE) - virgin

- Conforms to FDA CFR 177.1520
- BSE/TSE free

Closure Type: 2 steel wires (non product contact)

- Covered with FDA approved pressure sensitive tape with a backing of polypropylene or PVC

Bag Dimensions: 114 x 382mm

Thickness of Bag: 62 microns

Nominal Bag Volume: 1200ml

Method of Sterilisation: Sterility is assured due to extrusion temperatures exceeding 220°C which eliminates any present bacteria. During manufacturing all bags are sealed at both top and bottom in a continuous format before being cut into individual bags. This method assures that the bag's interior will never come in contact with its surroundings. Periodic validation tests, based on USP 34 – Sterility testing, are performed in the manufacturing facility to guaranty the sterility of the processes.

Packaging: 1000 bags per case (2 x 500)

Recommended Storage Conditions: Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture

ATEX: This product has not been ATEX tested

3. Country of Origin: Made and packed in Canada.**4. Repeated Use:** This is a single-use, disposable product and should not be re-used.**5. Approval:**

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 21-Jan-25

Note: The information given in this document has been given to us by our suppliers and is passed on in good faith, we believe the information to be correct.