


1. Part Number: 8175H-200S

2. Photograph:



3. Product Information:

Description:	SteriWare® Pallet Knife - 200mm – Sterile
Use:	Single use, disposable Pallet Knife
Material of Construction:	High Density Polyethylene (HDPE) - virgin • Conforms to FDA 21 CFR 177.1520 • Conforms to EU Regulations 10/2011 • Conforms to EC Regulations 1935/2004 • Conforms to REACH regulation 1907/2006 • BSE/TSE free
Method of Construction:	Single piece, injection moulded
Overall Length:	310mm
Blade Length:	200mm
Blade Width:	45mm
Nominal Weight:	55g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 50kGy)
Quantity per Box:	100
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This product has not been ATEX tested

4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

change.control@sampling.com

5. Country of Origin:

Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable product. Repeated use will lead to failure.

7. Approval:

Signed:

**Name:** Kuldeep Flora**Position:** Quality Manager**Date:** 30-Jul-24