

1. Part Numbers: 8180A-1000

2. Photograph:



3. Product Information:

Description	PP Tub 1000ml
Material of Construction:	Polypropylene (PP)
Nominal Volume:	1000ml
Nominal Height:	127mm
Max. Lid Diameter:	113mm
Internal Diameter:	104mm
Nominal Weight:	94g
Graduated:	No
Individually Bagged:	No
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 30-Jul-24