

1. Part Numbers: 8180A-25

2. Photograph:



3. Product Information:

Description	PP Tub 25ml
Material of Construction:	Polypropylene (PP)
Nominal Volume:	25ml
Brimful Volume:	25ml
Nominal Height:	50mm
Max. Lid Diameter:	30mm
Internal Diameter:	24mm
Nominal Weight:	7g
Graduated:	No
Individually Bagged:	No
BSE/TSE Free:	Yes
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 30-Jul-24