

**1. Part Number:** 8220A-003

**2. Label Design:**



The label is yellow with a white border. At the top, it says "Quality Control" in large white letters. Below that, it says "Contrôle Qualité" and "Control de Calidad" in smaller white letters. The main body of the label is divided into five horizontal sections by thin yellow lines. On the left side, there is vertical text: "Sampling Systems - Tel: +44 (0)1675 466992 E-mail: info@sampling.co.uk". On the right side, there is vertical text: "Sampling Systems Ltd. - www.sampling.co.uk - Reorder Code: 8220A-003".

**3. Product Information:**

|                                 |                                                         |
|---------------------------------|---------------------------------------------------------|
| Description:                    | Quality Control Label - Yellow - 150 x 150mm            |
| Adhesive:                       | Ultra Adhesive - High strength foil backed paper labels |
| Colours:                        | Yellow printing on a white label                        |
| FDA Acceptable Adhesive?        | No                                                      |
| Dimensions:                     | 150 x 150mm                                             |
| Number of Labels per Roll:      | 250 (approx).                                           |
| Weight:                         | 1300g per roll                                          |
| Method of Sterilisation:        | N/A                                                     |
| Recommended Storage Conditions: | Dry and ambient temperature                             |

**4. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Aug-24