



1. Part Number: 8220A-005

2. Label Design:

Quality Control	
Contrôle Qualité	Control de Calidad

Sampling Systems Ltd. - www.sampling.co.uk - Reorder Code 8220A-005

Sampling Systems - Tel: +44 (0)1675 466992 E-mail: info@sampling.co.uk

3. Product Information:

Description:	Quality Control Label - Black - 150 x 150mm
Adhesive:	Ultra Adhesive - High strength foil backed paper labels
Colours:	Black printing on a white label
FDA Acceptable Adhesive?	No
Dimensions:	150 x 150mm
Number of Labels per Roll:	250 (approx).
Weight:	1300g per roll
Method of Sterilisation:	N/A
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Aug-24