


1. Part Number: 8220A-32

2. Label Design:



The image shows a rectangular orange label with white text and input fields. At the top, the word "QUARANTINE" is written in large, bold, white capital letters. Below this, there are several fields for data entry: "Container" followed by two small white boxes and the word "of", "Material" followed by a long white box, "Batch / Lot Number" followed by a long white box, "Expiry Date" and "Part Number" each followed by a white box, and "Quantity" and "Initials & Date" each followed by a white box. On the left side of the label, there is vertical text: "Sampling Systems Ltd. Tel +44 (0)1675 466992 E-mail: info@sampling.co.uk". On the right side, there is vertical text: "Re-order code 8220A-32 Visit our website at: www.sampling.co.uk".

3. Product Information:

Description:	'Quarantine' Label - 95 x 95mm
Adhesive:	Ultra Adhesive - High strength foil backed paper labels
Colours:	Orange printing on a white label
FDA Acceptable Adhesive?	No
Dimensions:	95 x 95mm
Number of Labels per Roll:	500 (approx).
Weight:	1010g per roll
Method of Sterilisation:	N/A
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 6-Aug-24