

1. Part Number: 8333P-2

2. Photograph:



3. Product Information:

Description: Bailer 2000ml
Material of Construction: Polypropylene (PP)

- FDA 21 CFR 177.1520
- EU 10/2011
- EC 1935/2004
- BSE/TSE free

Construction: Hygienic, single piece construction
Nominal Volume: 2000ml
Overall Length: 355mm
Handle Length: 135mm
Bowl Diameter: 220mm
Bowl Height: 79mm
Nominal Weight: 280g
Recommended Storage Conditions: Dry and ambient temperature

Note:

This device is manufactured in a general factory environment. It is designed to be cleaned and re-used. For hygienic applications this product must be washed before first use.

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Aug-24