

**1. Part Number:** 8370B-100S

**2. Photograph:**



### 3. Product Information:

Description:	SteriWare Funnel 100mm diameter - PP Blue - Sterile
Use:	Single use, disposable liquid funnel
Material of Construction:	Polypropylene (PP) - virgin <ul style="list-style-type: none"> <li>• Conforms to FDA 21 CFR 177.1520</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• Conforms to REACH regulation 1907/2006</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Max. Diameter:	117mm
ID at Top of Funnel:	100mm
Overall Height:	103mm
OD at Bottom of Spout:	11mm
ID at Bottom of Spout:	8mm
OD at Top of Spout:	17mm
Nominal Weight:	17g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25 to 45 kGy)
Number of Samplers per Box:	50
Recommended Storage Conditions:	Dry and ambient temperature

Shelf Life: 5 years from the date of manufacture  
ATEX: This product has not been ATEX tested

**4. Change Control:** Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:** Made and packed in the UK.

**6. Repeated Use:** This is a single-use, disposable sampler. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Aug-24