

**1. Part Number:** 8404H-250S

**2. Photograph:**



### 3. Product Information:

Description	Wide Mouth HDPE Bottle 250ml - Sterile
Materials of Construction:	Bottle HDPE (conforms to FDA 177.1520 & USP Class VI) Cap PP (conforms to FDA 177.1520 & USP Class VI)
Manufacturing Environment:	Controlled environment within factory
Packing Environment:	Class 10,000 Cleanroom
Nominal Volume:	250ml
Brimful Volume:	290ml
Nominal Height:	131mm
Diameter:	62mm
Internal Neck Diameter:	32mm
Weight:	37g
Individually Bagged:	Yes (in a PE/PET laminate film bag)
Quantity per box:	50
Method of Sterilisation:	Gamma Irradiation (25-45 kGy)
BSE/TSE Free:	Yes
Supplied with lid screwed on:	Yes
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture
ATEX:	This device has not been ATEX tested

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Packed and sterilised in the UK.

**6. Repeated Use:**

This is a single-use, disposable device. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 22-Aug-24