

**1. Part Number:** 8703A-05

**2. Photograph:**



### 3. Product Information:

Description: Titrating Automatic Burette 5ml

Product Contact Materials:

- Polyethylene
- Borosilicate Glass
- Latex

Method of Construction: Parts push together, they can be dis-assembled for cleaning

Overall Height: 480mm

Capacity of Burette: 5ml

Graduation Marks on Burette: Every 0.05ml  
(Accuracy +/- 0.03ml, Category B ISO 385)

Volume of Reservoir Bottle: 500ml

Nominal Weight: 147g

ATEX: This device has not been ATEX tested

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25