

1. Part Number: 8703A-15

2. Photograph:



3. Product Information:

Description:	Titrating Automatic Burette 15ml
Product Contact Materials:	<ul style="list-style-type: none">- Polyethylene- Borosilicate Glass- Latex
Method of Construction:	Parts push together, they can be dis-assembled for cleaning
Overall Height:	490mm
Capacity of Burette:	15ml
Graduation Marks on Burette:	Every 0.10ml (Accuracy +/- 0.045ml, Category B ISO 385)
Volume of Reservoir Bottle:	500ml
Nominal Weight:	146g
ATEX:	This device has not been ATEX tested

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25