

1. Part Number: 8703A-25

2. Photograph:



3. Product Information:

| | |
|------------------------------|---|
| Description: | Titrating Automatic Burette 25ml |
| Product Contact Materials: | <ul style="list-style-type: none"> - Polyethylene - Borosilicate Glass - Latex |
| Method of Construction: | Parts push together, they can be dis-assembled for cleaning |
| Overall Height: | 610mm |
| Capacity of Burette: | 25ml |
| Graduation Marks on Burette: | Every 0.10ml (Accuracy +/- 0.045ml, Category B ISO 385) |
| Volume of Reservoir Bottle: | 1000ml |
| Nominal Weight: | 194g |
| ATEX: | This device has not been ATEX tested |

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 9-Jan-25