

**1. Part Number:** 8703A-50

**2. Photograph:**



### 3. Product Information:

Description:	Titrating Automatic Burette 50ml
Product Contact Materials:	<ul style="list-style-type: none"> <li>- Polyethylene</li> <li>- Borosilicate Glass</li> <li>- Latex</li> </ul>
Method of Construction:	Parts push together, they can be dis-assembled for cleaning
Overall Height:	870mm
Capacity of Burette:	50ml
Graduation Marks on Burette:	Every 0.10ml (Accuracy +/- 0.075ml, Category B ISO 385)
Volume of Reservoir Bottle:	1000ml
Nominal Weight:	225g
ATEX:	This device has not been ATEX tested

### 4. Approval:

Signed:



Name: Kuldeep Flora  
Position: Quality Manager  
Date: 9-Jan-25