

**1. Part Number:** 8890P-250S

**2. Photograph:**



**3. Product Information:**

Description:	SteriWare® Powder Bottle 250ml - Sterile
Use:	Single use, disposable bottle
Material of Construction:	Polypropylene (PP) – opaque <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1520</li> <li>• Conforms to EU Regulation 10/2011</li> <li>• Conforms to EC Regulation 1935/2004</li> <li>• BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Nominal Volume:	250ml
Brimful Volume:	280ml
Method of Construction:	Single piece, injection moulded
Height:	88mm
Lid Diameter:	76mm
Nominal Weight:	56g
Moulding Environment:	Class 100,000 Cleanroom
Primary Packing Environment:	Class 10,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)

Method of Sterilisation: Gamma Irradiation (25-45kGy)  
Number of Bottles per Box: 50  
Recommended Storage Conditions: Dry and ambient temperature  
Shelf Life: 5 years from the date of manufacture  
ATEX: This product has not been ATEX tested.

#### 4. Change Control:

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to

[change.control@sampling.com](mailto:change.control@sampling.com)

#### 5. Country of Origin:

Made and packed in the UK.

#### 6. Repeated Use:

This is a single-use, disposable device. Repeated use will lead to failure.

#### 7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-24