

**1. Part Number:** 8916A-100

**2. Photograph:**



### 3. Product Information:

|                                 |                                                                                                                                                                                                                                                      |
|---------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Description:                    | SteriWare® Scoop with Lid 100ml                                                                                                                                                                                                                      |
| Use:                            | Single use, disposable scoop with transparent click on lid                                                                                                                                                                                           |
| Nominal Flat-filled Volume:     | 100ml                                                                                                                                                                                                                                                |
| Material of Construction:       | High Impact Polystyrene (HIPS) – rigid, white <ul style="list-style-type: none"> <li>• Conforms to FDA CFR 177.1640</li> <li>• Conforms to EU Regulations 10/2011</li> <li>• Conforms to EC Regulations 1935/2004</li> <li>• BSE/TSE free</li> </ul> |
| Method of Construction:         | Single piece, injection moulded                                                                                                                                                                                                                      |
| Scoop Overall Length            | 200mm                                                                                                                                                                                                                                                |
| Scoop Handle Length:            | 95mm                                                                                                                                                                                                                                                 |
| Scoop Maximum Width:            | 59mm                                                                                                                                                                                                                                                 |
| Nominal Weight:                 | 37g (Scoop 25g + Lid 12g)                                                                                                                                                                                                                            |
| Moulding Environment:           | Class 100,000 Cleanroom                                                                                                                                                                                                                              |
| Primary Packing Environment:    | Class 10,000 Cleanroom                                                                                                                                                                                                                               |
| Individually Bagged?            | Yes (heat sealed PE/PET bag)                                                                                                                                                                                                                         |
| Method of Sterilisation:        | N/A                                                                                                                                                                                                                                                  |
| Number of Samplers per Box:     | 50                                                                                                                                                                                                                                                   |
| BSE/TSE Free:                   | Yes                                                                                                                                                                                                                                                  |
| Recommended Storage Conditions: | Dry and ambient temperature                                                                                                                                                                                                                          |
| Shelf Life:                     | 5 years from the date of manufacture                                                                                                                                                                                                                 |
| ATEX:                           | This device has not been ATEX tested                                                                                                                                                                                                                 |

**4. Change Control:**

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to [change.control@sampling.com](mailto:change.control@sampling.com)

**5. Country of Origin:**

Made and packed in the UK.

**6. Repeated Use:**

This is a single-use, disposable scoop. Repeated use will lead to failure.

**7. Approval:**

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 3-Sep-24