

1. Part Number: A104-300

2. Photograph:



3. Product Information:

Description:	Double Ended Spoon 300mm
Material of Construction:	Stainless steel
Overall Length:	300mm
Bowl Sizes:	2.0ml & 8.0ml
Nominal Weight:	57g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 16-Dec-24