




1. Part Number: A126-200

2. Photograph:



3. Product Information:

Description

Stainless Storage Drum 200L

Materials of Construction:

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Handles: 304 stainless steel
- Gasket: Silicone (FDA 177.2600 compliant & BSE/TSE free)

Nominal Volume:

200 Litres

Brimful Volume:

204 Litres

Method of Construction:

Welds ground and polished

Surface Finish:

Better than 0.5 microns Ra

Dimensions

A = 570 mm
B = 590 mm
C = 800 mm

Drum Wall Thickness:

1.2mm

Graduation Marks:

None

Handles:

2 side handles welded onto patch plates

Weight:

20.2 Kg (body + gasket + lid + clamp)

BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Recommended Storage Conditions:

Dry and ambient temperature

Notes:

These containers are primarily designed for storing powders. They are not suitable for the storage or transportation of beer.

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 20-Sep-24