

1. Part Number: A236-350

2. Photograph:



3. Product Information:

Description	Deep Tray
Nominal Volume:	4L
Material of Construction:	316L stainless steel
Method of Construction:	Single piece, crevice free
Surface Finish:	Better than 0.5 microns Ra
Length:	380mm
Width:	280mm
Height:	54mm
Nominal Weight:	905g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24