

1. Part Number: A296-1

2. Photograph:



3. Product Information:

| | |
|---------------------------------|--|
| Description | Stainless Steel Bailer 1L |
| Materials of Construction: | 316L stainless steel |
| Nominal Volume: | 1 Litre |
| Method of Construction: | Welds ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Handle Length: | 115mm |
| Max. Bowl Diameter: | 170mm |
| Bowl Height: | 70mm |
| Nominal Weight: | 435g |
| Graduation Marks: | None |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin. |
| Recommended Storage Conditions: | Dry and ambient temperature |

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 18-Dec-24