

1. Part Number: A316-5000

2. Photo:



3. Product Information:

| | |
|---------------------------------|---|
| Description | PharmaScoop® 5000ml |
| Nominal Volume: | 5000ml |
| Material of Construction: | 316L stainless steel |
| Method of Construction: | Crevice free, all welds ground and polished |
| Surface Finish: | Better than 0.5 microns Ra |
| Overall Length: | 355mm |
| Body Length: | 225mm |
| Body Diameter: | 220mm |
| Nominal Weight: | 1165g |
| Total Surface Area: | 2.55x10 ⁵ mm ² (approx.) |
| BSE/TSE statement: | All polishing compounds used in the manufacture of this product are of vegetable origin |
| Recommended Storage Conditions: | Dry and ambient temperature |

Note: This product should be thoroughly cleaned before it is used for the first time

4. Approval:

Signed:

Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24