

**1. Part Number:** A336-4000

**2. Photograph:**



### 3. Product Information:

Description:	Open PharmaScoop® 1500ml
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	4000ml
Overall Length:	430mm
Body Length:	310mm
Body Diameter:	205mm
Scoop Height:	85mm
Nominal Weight:	950g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.