

1. Part Number: A336-750**2. Photograph:****3. Product Information:**

Description:	Open PharmaScoop® 750ml
Material of Construction:	316L Stainless Steel
Method of Construction:	Crevice free, all interior welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Nominal Volume:	750ml
Overall Length:	320mm
Body Length:	205mm
Body Diameter:	125mm
Scoop Height:	50mm
Nominal Weight:	420g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 18-Dec-24

Note: All dimensions are nominal. Please contact us if any dimension is critical for your application.