


**1. Part Number:** A566-02

**2. Photograph:**



A = 119mm

B = 98mm

C = 185mm

D = 135mm

E = 130mm

**3. Product Information:**

**Description:**

TC Container 2L

**Materials of Construction:**

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)

**Method of Construction:**

All welds ground and polished  
Crevice free interior

**Surface Finish:**

Better than 0.5 microns Ra

**Nominal Volume**

25L

**BSE/TSE statement:**

All polishing compounds used in the manufacture of this product are of vegetable origin.

**Handles:**

Two side handles welded on to patch plates

**Graduation Marks:**

None

**Pressure Rating:**

None

**Weight:**

2.1kg (body + lid + gasket + clamp)

**Recommended Storage Conditions:**

Dry and ambient temperature

**Notes:**

The neck of the container is fitted with a 4" diameter ferrule which conforms to BS4825-3

These containers are primarily designed for storing powders. They are not suitable for the storage or transportation of beer.

## 4. Spare Gaskets:

Supplied with a 4" silicone gasket as standard. Gaskets are consumable parts and will need to be inspected and changed at the first signs or wear.

Part No.	Description
K201-4	Gasket 4" silicone – FDA acceptable grade
K202-4	Gasket 4" EPDM – FDA acceptable grade
K203-4	Gasket 4" Viton – FDA acceptable grade

## 5. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 13-Jan-26