

1. Part Number: A606-750

2. Photograph:



3. Product Information:

Description	GMP PharmaScoop® 750ml
Nominal Volume:	750ml
Material of Construction:	316L stainless steel
Method of Construction:	Crevice free, all welds ground and polished
Surface Finish:	Better than 0.5 microns Ra
Overall Length:	333mm
Body Length:	180mm
Body Diameter:	98mm
Nominal Weight:	465g
BSE/TSE statement:	All polishing compounds used in the manufacture of this scoop are of vegetable origin
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 7-Jan-25