

1. Part Number: A625-280

2. Photograph:



3. Product Information:

Description:	Chemical Spoon 200mm
Material of Construction:	304 Stainless steel
Overall Length:	280mm
Bowl Length:	90mm
Bowl Width:	50mm
Nominal Bowl Volume:	25ml
Thickness of metal:	1.5mm
Nominal Weight:	108g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25