

1. Part Number: A625-320

2. Photograph:




3. Product Information:

Description:	Chemical Spoon 320mm
Material of Construction:	304 Stainless steel
Overall Length:	317mm
Bowl Length:	110mm
Bowl Width:	48mm
Nominal Bowl Volume:	35ml
Thickness of metal:	1.5mm
Nominal Weight:	126g
BSE/TSE statement:	All polishing compounds used in the manufacture of this product are of vegetable origin.
Recommended Storage Conditions:	Dry and ambient temperature

4. Approval:

Signed:



Name: Kuldeep Flora
Position: Quality Manager
Date: 7-Jan-25