


**1. Part Number:** A746-03

**2. Photograph:**



### 3. Product Information:

**Description:**

Jacketed Container 3L

**Materials of Construction:**

- Body: 316L stainless steel
- Lid: 316L stainless steel
- Clamp: 304 stainless steel
- Gasket: Silicone (FDA acceptable grade)

**Method of Construction:**

All welds ground and polished  
Crevice free interior

**Surface Finish:**

Better than 0.5 microns Ra

**Nominal Volume:**

3 Litres

**Dimensions:**


A = 207mm  
B = 160mm  
C = 190mm  
D = 40mm  
E = 115mm  
F = 1.2mm  
G = 119mm

**Neck / Lid Size:**

4" BS4825-3

**Jacket Inlet/Outlet Size:**

Compatible with 1" and 1.5" tri-ferrule connections BS4825-3



BSE/TSE statement:

All polishing compounds used in the manufacture of this product are of vegetable origin.

Weight:

3.4 Kg (body + lid + gasket + clamp)

Recommended Storage Conditions:

Dry and ambient temperature

#### 4. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 8-Jan-25